

INSUFLATION FILTER

FLEXIBLE TRANSPARENT TUBE

Flexible and lightweight tubing helps prevent folds

LUER LOCK CONNECTOR

BIDIRECTIONAL HYDROPHOBIC FILTERING MEMBRANE

of 0.2 μ m and 99.99%
filtration efficiency

SUCTION CONNECTOR

Compatible with all standard insufflators

Filter with a 2.5m flexible tube designed for CO₂ filtration with 99.99% efficiency, preventing cross-contamination between patient and insufflator.



WWW.AMERICANINSTRUMENTS.US

CONTACT: +1 (407)954-0004



American Instruments

