

LAPAROHEPA FILTER

LUER LOCK CONNECTOR

**ULPA FILTER +
ACTIVATED CARBON MEMBRANE**

Removes 99.999% of 0.1
micron particles +
activated carbon membrane
that absorbs odors
and chemical toxins.

TRANSPARENT FLEXIBLE TUBE

The LaparoHepa prevents contamination of the surgical environment during laparoscopic procedures. It eliminates biological pathogens (bacteria and viruses) as well as chemical toxins and particles resulting from cauterization.



WWW.AMERICANINSTRUMENTS.US

CONTACT: +1 (407)954-0004



American Instruments

