

# MAGNETIC MAT



Magnetic Mat helps prevent items from falling.

Designed with strong round-raised magnets, the silicone Magnetic Mat measurements are 30 x 40 cm and helps prevent items from falling, allowing healthcare professionals to grasp for easy removal.

Intended to be utilized as a hands-free “neutral-zone” sharps transfer area, the reusable waterproof mat conforms to all surface shapes and can be fastened to drapes with a standard towel clamp for ease of access.



REF: AI0111 - MANTA MAGNÉTICA



American Instruments