

# SAFE BAG<sup>®</sup>

**ACCESS 3**  
Morcellator

**ACCESS 2**  
Auxiliary instruments

**ACCESS 1**  
Optic

**Measuring marks**  
Graphic marks to guide surgeons  
in measuring tissues sizes and bulk

**Material**  
Made by translucent and  
high resistant material

SAFE BAG is a medical device used in minimally invasive surgery that allows an efficient protection barrier against cell dispersion into abdominal cavity.

This medical device should only be used by trained and licensed professional for minimally invasive surgery (videolaparoscopy).



[WWW.AMERICANINSTRUMENTS.US](http://WWW.AMERICANINSTRUMENTS.US)

CONTACT: +1 (407) 954-0004



American Instruments

